


Agrimin 24-7 Copper Boluses for Cattle





Bolus Technology

Copper oxide rods held together in wax matrix

Wax dissolves in rumen and copper oxide rods released

Copper oxide rods slowly dissolve in abomasum

Active Component

Each Bolus Contains

Copper - 20g bolus

(as copper oxide rods)

Copper - 30g bolus (as copper oxide rods)

Copper - 36g bolus (as copper oxide rods)

Benefits

(20g copper oxide)

25.2g

16.8q

(30g copper oxide)

30.2g

(36g copper oxide)

• 0


- Optimise growth rates in youngstock and finishing animals
- Optimise coat condition
- Easy to administer
- No capsule means no breakages or damage prior to or during administration

Administration

Administer orally.
Place the bolus at the back of the mouth using a recommended
Agrimin bolus applicator:

- 20g bolus one bolus per 200kg liveweight
- 30g bolus one bolus per 300kg liveweight
- 36g bolus one bolus per 360kg liveweight
- Available in packs of 24 boluses

When to use Agrimin 24-7 Copper





35 Beach Street, Kippa Ring, QLD 4021

Tel: +61 7 3283 5077 Fax: +61 7 3283 5088 Email: rural@pacificbiologics.com.au Web: www.pacificbiologics.com.au











AUSA4 - 16/2